

Technical Information

Autosampler Compatibility Table

This table indicates the categories of vials that are compatible with various models of autosamplers. Certain autosamplers require the purchase of optional vial trays and, in few cases, programming upgrades to use all of the vials listed.

An X mark indicates that the vials from this category are compatible with the autosampler in most configurations.

An F mark indicated that a cap having an outer flange is required for the vial to operate properly with the autosampler.

An M mark indicates that a magnetic seal is required for use with the autosampler.

MANUFACTURER	MODEL	8mm CRIMP	11mm CRIMP	8-425 SCREW THREAD	TARGET DP SCREW THREAD	TARGET 10-425 SCREW THREAD	15X45mm 4mL	HEADSPACE	WEBSEAL
Agilent Technologies®	1050	x	x		x				
	1100/1200		x		x				
	1090A	x	x	x	x				
	1100/1200 Wellplate/ NanoFlow		x		x				x
	6850ALS		x						
	7673A/7683		x						
	7694E							x	
	220								x
AI	42 Position Tray		x						
	60 Position Tray	x							
AIM	CPS-100/200		x	x	x				
Alcott®	708AL								
	718AL		x	x	x				x
Alpha M.O.S.	Fox/Prometheus		x	x	x			M	
	Kronos		x	x	x				
Antec Leyden	ALEXYS AS 100		x	x	x				
ATAS GL	Focus		x	x	x			M	
Beckman Coulter®	501/507/System Gold	x	x	x	x		x		
Bruker	Map II		x	x	x	x			x
	LC51						x		
CTC™	see Leap Technologies		x	x	x				x
DANI	ALS-1000/39.80/86.80		x						
	HSS 39.50/ 86.50/ SPT 37.50							x	
Dionex	ASI-100/A550		x	x	x	x	x		
	AS-40			x			x		
	FAMOS		x	x	x				x
DynaTech	LC-241	x	x	x	x				
ESA	540		x	x	x	x	x		
	540-MT		x	x	x	x			x
Fisons®	A-200S/AS 150/800/8000		x	x	x				
	AS 200	x							
GBC	LC 1650		x	x	x	x	x		
Gerstel	MPS 2	x	x	x	x	x		M	x
	MPS 3	x	x	x	x	x			
Gilson®	234		x	x	x				
	235		x	x	x				x
	231 XL/232 XL/233 XL	x	x	x	x				
	215		x	x	x	x			x
	250		x	x	x				x
Jasco®	851/AS950/1550/1555			F					
	AS-2059/AS2059Plus	x		F		x			x
	AS-2055/AS/2057	x		F		x			x
LEAP Technologies	CLPAL/GC PAL	x	x	x	x	x			
	Combi PAL	x	x	x	x	x			
	HTC PAL/HTS PAL	x	x	x	x	x			

MANUFACTURER	MODEL	8mm CRIMP	11mm CRIMP	8-425 SCREW THREAD	TARGET DP SCREW THREAD	TARGET 10-425 SCREW THREAD	15X45mm 4mL	HEADSPACE	WEBSEAL
Hitachi®	AS-1000		x	x	x				
	AS-2000		x	x	x				
	AS-4000		x		x				
	Lachrom L-7200/7250		x	x	x				
	Lachrom Elite L-2200		x	x	x				x
Metrohm	Triathlon		x	x	x	x			
Perkin Elmer™	Autosystem/AS-2000/XL		x	x	x		x		
	Clarus 500/600		x		x				
	Integral 4000		x	x	x				
	TurboMatrix 40/110							x	
	ISS-100/200/Series 200	x	x	x	x				
	HS16/40							x	
	LC Plus	x	x	x	x				x
Polymer Labs	GPC 110/210/220		x						
Selerity	3100		x	x	x	x			
Shimadzu®	AOC-8B/9, AOC-14/1400		x	x	x				
	AOC-5000	x	x	x	x	x		M	x
	AOC-20		x		x		x		
	HSS-2B/4B							x	
	SIL-6A/6B/9A			F		x			
	SIL-10A			x	x	x			
	SIL-HT/10ADVP	x	x	x	x	x			x
LC-2010	x	x	x	x	x			x	
Spark-Holland	Symbiosis/Reliance		x	x	x	x			x
Thermo Scientific	Surveyor LC		x	x	x				x
	AS 3000/Trace GC/Focus LC		x		x				
	TriPlus		x	x	x	x		x	
Tosoh®	TSK-6080/AS-8010/AS-8020		x	F		x	x		
Teledyne/Tekmar®	7000/7050							x	
	HT3A							x	
Varian®	CP-8410		x	x	x	x	x		
	CombiPAL	x	x	x	x	x		M	
	8034/8035/8100/8200		x	x	x				
	9095/9100		x	x	x				
	Prostar 400/410/420		x	x	x	x			
	Prostar 430		x	x	x	x			x
	GeneSis							x	
Waters™	WISP 96 Pos. Tray	x							
	WISP 48 Pos. Tray						x		
	717Plus	x	x						
	Alliance 2690			x	x	x			
	Model 2767			x	x	x			x
	Model 2777					x			x
	Acquity		x	x	x	x			x
CapLC		x	x	x				x	